

AVENIO ctDNA and Tumor Tissue Expanded Panel



This 77-gene panel, which comes with the AVENIO ctDNA Expanded Kit and the AVENIO Tumor Tissue Expanded Kit, contains 17 biomarkers in the U.S. National Comprehensive Cancer Network (NCCN) and other guidelines, plus 60 additional biomarkers investigated in clinical trials.*

Gene	All Coding Regions	CRC	Lung	Breast	Gastric	Prostate	Bladder	Glioma	Melanoma	Ovarian	Thyroid	Pancreatic	Other Cancers	SNV	Indel [‡]	Fusion [§]	CNV [§]
ABL1													CML, ALL	●			
AKT1		●	●	●						●				●			
AKT2		●	●	●				●		●		●	Thymoma	●			
ALK			●										ALCL, Neuroblastoma, Inflammatory Myofibroblastic Tumor, Spitzoid Tumor, Large B-Cell Lymphoma, Rhabdomyosarcoma	●	●	●	
APC		●			●			●				●	HCC, Medulloblastoma, Desmoid	●	●		
AR	●			●		●								●			
ARAF			●											●			
BRAF		●	●					●	●	●	●		Borderline Ovarian, Cholangiocarcinoma, Spitzoid Tumor, Pancreas Acinar Carcinoma, Melanocytic Nevus, GIST	●	●		
BRCA1	●			●						●		●		●			
BRCA2	●			●				●		●		●		●			
CCND1	●	●		●									CLL, B-ALL, MM, Mantle Cell Lymphoma	●			
CCND2	●												NHL, CLL	●			
CCND3	●						●	●					MM, Burkitt Lymphoma, other B-Cell Lymphomas	●			
CD274	●												Primary Mediastinal B-Cell Lymphoma, Hodgkin Lymphoma, CML, RCC, Lymphoepithelial-like Carcinoma	●			
CDK4	●			●					●				Liposarcoma	●			
CDK6				●				●					Endometrial Squamous Cancer, T-Cell Leukemia	●			
CDKN2A	●								●			●	Astrocytoma	●			
CSF1R				●		●				●			Endometrial, Testicular GCT, Malignant Histiocytosis	●			
CTNNB1		●								●			HCC, Medulloblastoma, Mesothelioma	●	●		

All coding regions are based on the longest transcript from Ensembl build 82.

* National Comprehensive Cancer Network. <http://www.nccn.org>. October 15, 2016.

‡ Indels are limited to variants in a pre-specified list of positions, referred to as "Loci of Interest", except for EGFR exon 19 long deletions, EGFR exon 20 long insertions and MET long insertions, which are not restricted to a pre-defined set of Indels.

§ Detection of Fusions and CNVs are limited to variants in a pre-specified list of positions, referred to as "Loci of Interest" in the AVENIO analysis software.

Genes in the AVENIO ctDNA and Tumor Tissue Expanded Panel



Gene	All Coding Regions	CRC	Lung	Breast	Gastric	Prostate	Bladder	Glioma	Melanoma	Ovarian	Thyroid	Pancreatic	Other Cancers	SNV	Indel [‡]	Fusion [§]	CNV [§]
DDR2			•										Lung Adenoid Cystic Carcinoma	•			
DPYD**		•												•			
EGFR	•		•					•						•	•		•
ERBB2	•	• [†]	•	•	•			•		•				•	•		•
ESR1	•			•										•			
EZH2													MDS, Mantle Cell Lymphoma, CML, DLBCL, Cholangiolocellular Carcinoma	•			
FBXW7	•	•		•									Endometrial, T-ALL, Lymphoblastic Leukemia	•			
FGFR1			•	•									MPN, NHL, Salivary Adenoma	•			
FGFR2			•	•	•								Endometrial	•		•	
FGFR3			•				•						MM, T-Cell Lymphoma, Cervical, Testicular GCT	•		•	
FLT1													CNS Hemangiomas	•			
FLT3													AML, ALL	•			
FLT4													Soft Tissue Sarcoma	•			
GATA3				•									RCC	•			
GNA11									•					•			
GNAQ									•					•			
GNAS		•											Pituitary Adenoma	•			
IDH1								•					AML	•			
IDH2								•					AML	•			
JAK2													ALL, AML, MPN, CML	•			
JAK3													Acute Megakaryocytic Leukemia, ETP ALL	•			
KDR			•										Angiosarcoma	•			
KEAP1	•		•											•			
KIT									•				AML, Testicular Germ Cell Tumors, GIST, Thymic Carcinoma	•	•		

All coding regions are based on the longest transcript from Ensembl build 82.

** Inactivating mutations are associated with sensitivity to fluoropyrimidines.

[†] Kavuri et al. Cancer Discovery 8: 832-41 (2015). PMID: 26243863.

[‡] Indels are limited to variants in a pre-specified list of positions, referred to as "Loci of Interest", except for EGFR exon 19 long deletions, EGFR exon 20 long insertions and MET long insertions, which are not restricted to a pre-defined set of Indels.

[§] Detection of Fusions and CNVs are limited to variants in a pre-specified list of positions, referred to as "Loci of Interest" in the AVENIO analysis software.

Genes in the AVENIO ctDNA and Tumor Tissue Expanded Panel



Gene	All Coding Regions	CRC	Lung	Breast	Gastric	Prostate	Bladder	Glioma	Melanoma	Ovarian	Thyroid	Pancreatic	Other Cancers	SNV	Indel [‡]	Fusion [§]	CNV [§]
KRAS	●	●	●		●					●	●	●	AML, Juvenile Myelomonocytic Leukemia	●			
MAP2K1		●	●						●					●			
MAP2K2			●						●					●			
MET	●		●										Head and Neck Squamous Cell, HCC, Papillary RCC	●	●		●
MLH1	●	●		●	●					●			Endometrial, Small Intestine	●			
MSH2	●	●		●	●					●			Endometrial, Small Intestine	●			
MSH6	●	●		●	●					●			Endometrial, Small Intestine	●			
MTOR								●					Endometrial, Head and Neck SCC, Clear Cell RCC	●			
NF2	●												Meningioma, Acoustic Neuroma, Mesothelioma	●			
NFE2L2			●										Head and Neck SCC	●			
NRAS		●	●						●		●		MM, AML, Juvenile Myelomonocytic Leukemia	●			
NTRK1			●								●		Spitzoid Tumor	●		●	
PDCD1LG2	●												PMBL, Hodgkin Lymphoma	●			
PDGFRA													GIST, Hypereosinophilic Syndrome	●			
PDGFRB													MPN, AML, CMML, CML, Juvenile Myelomonocytic Leukemia	●			
PIK3CA		●		●	●			●		●	●		HCC	●	●		
PIK3R1		●						●		●				●			
PMS2	●	●		●	●					●			Endometrial, Small Intestine	●			
PTCH1													Basal Cell Carcinoma, Medulloblastoma	●			
PTEN	●	●		●		●		●			●		Endometrial, Head and Neck SCC, RCC	●	●		
RAF1						●							Pilocytic Astrocytoma	●			
RB1	●			●			●						Retinoblastoma, Small Cell Lung Carcinoma, Osteosarcoma	●			

All coding regions are based on the longest transcript from Ensembl build 82.

[‡] Indels are limited to variants in a pre-specified list of positions, referred to as "Loci of Interest", except for EGFR exon 19 long deletions, EGFR exon 20 long insertions and MET long insertions, which are not restricted to a pre-defined set of Indels.

[§] Detection of Fusions and CNVs are limited to variants in a pre-specified list of positions, referred to as "Loci of Interest" in the AVENIO analysis software.

Genes in the AVENIO ctDNA and Tumor Tissue Expanded Panel



Gene	All Coding Regions	CRC	Lung	Breast	Gastric	Prostate	Bladder	Glioma	Melanoma	Ovarian	Thyroid	Pancreatic	Other Cancers	SNV	Indel [‡]	Fusion [§]	CNV [§]
RET		●	●								●		Pheochromocytoma, Spitzoid Tumor, Mucosal Neuroma	●		●	
RNF43		●								●		●	Cholangiocarcinoma	●			
ROS1			●		●			●					Spitzoid Tumor, Cholangiocarcinoma, Borderline Ovarian	●		●	
SMAD4	●	●										●	Small Intestine	●			
SMO	●												Basal Cell Carcinoma, Medulloblastoma	●			
STK11	●		●						●			●	Testicular GCT	●			
TERT promoter							●	●	●				HCC, Basal Cell Carcinoma, Skin Squamous Cell, Mesothelioma, Medulloblastoma	●			
TP53	●	●	●	●	●			●		●		●	Various Leukemias, Adenocortical, Basal Cell Carcinoma, Esophageal, Soft Tissue and Bone Sarcomas, Head and Neck SCC	●	● ^{§§}		
TSC1							●						RCC, Angiomyolipoma	●	●		
TSC2													RCC, Angiomyolipoma, Head and Neck SCC	●			
UGT1A1 ^{††}		●												●			
VHL	●												Hemangioblastoma, Pheochromocytoma, RCC	●			

All coding regions are based on the longest transcript from Ensembl build 82.
^{††} Certain variants are associated with sensitivity to irinotecan.
[‡] Indels are limited to variants in a pre-specified list of positions, referred to as "Loci of Interest", except for EGFR exon 19 long deletions, EGFR exon 20 long insertions and MET long insertions, which are not restricted to a pre-defined set of Indels.
[§] Detection of Fusions and CNVs are limited to variants in a pre-specified list of positions, referred to as "Loci of Interest" in the AVENIO analysis software.
^{§§} Indel for TP53 is called for tissue only, not plasma.

Associations of variations or mutations in these genes with the listed diseases were sourced from the following: COSMIC, GeneCards, My Cancer Genome. Analysis performed 11/2016.

台灣羅氏醫療診斷設備股份有限公司
104台北市中山區民權東路三段2號10樓
(02) 2183 6688
MC-TW-01069
僅供研究用 · 不可用於臨床診斷

